



12^{3/4}

29^{1/32}

29^{1/32}

9

3^{1/32}

3^{1/32}

PANTONE 326
PANTONE 425
PANTONE 1945

1/2

3

1

3

8^{15/32}

3^{1/32}

You're on your way to cleaner water!

Inhibits the growth of microbial contamination in dental waterlines.
Meets or exceeds <200 CFU/mL for dental unit waterlines.

IMS-1440

Active Ingredients:
Iodine*46.0%
Other Ingredients54.0%
TOTAL100.0%

This product contains 0.5 lb of iodine (I) per pound.

KEEP OUT OF REACH OF CHILDREN
CAUTION

*Bound to polystyrene divinyl benzene
quaternary ammonium anion exchange resin.



For Use as an Antimicrobial Agent in Dental Unit Waterlines
HU-FRIEDY WATERLINE FILTER
Manufactured for:
Hu-Friedy Mfg. Co., LLC
3232 N. Rockwell St.
Chicago, IL 60618
Net Contents: 1 Dental Waterline Filter
EPA Est. No. 73593-0H-2
EPA Reg. No. 75372-1
+D6715481113



HU-FRIEDY
WATERLINE FILTER
For Use as an Antimicrobial Agent in Dental Unit Waterlines

- Cleans and maintains dental unit waterlines
- Inhibits the growth of microbial contamination in dental waterlines
- Meets or exceeds <200 CFU/mL for dental unit waterlines
- Container has no silver or other hazardous waste – disposal as easy as throwing it in the trash
- Ability to be used with either municipal or independent reservoir
- One year life with minimal human interaction
- No daily, weekly or monthly protocols
- No corrosion or etching of metal parts
- No effect on bonding agents



FIRST AID

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

If on skin or clothing: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

If swallowed: Call a poison control center or doctor immediately for treatment advice.

Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give anything by mouth to an unconscious person.

Have the product container or label with you when calling poison control center or doctor, or going for treatment.

See sheet for additional precautionary statements, first aid and directions for use.

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS
CAUTION Causes moderate eye irritation. Avoid contact with eyes, skin and clothing. Wear protective eyewear (goggles) and gloves when handling. Wash thoroughly with soap and water after handling and before eating, drinking, gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

DIRECTIONS FOR USE

The exposed iodine resin pellets are toxic to fish and aquatic organisms. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge.

It is a violation of federal law to use this product in any manner inconsistent with its labeling, dental unit waterlines. The product reduces microbial contamination and odor causing bacteria in Hu-Friedy Waterline Filter is designed for use on dental unit waterlines attached to dynamic dental unit waterlines.

Prior to installing Hu-Friedy Waterline Filter, thoroughly clean dental unit waterlines following the manufacturer's instructions or with an appropriate cleaning agent. Insert Hu-Friedy Waterline Filter into the inlet and outlet fittings of dental unit waterlines.

Once the product is in place, the supply water will pass through the iodinated resin bed, which will continuously release low levels of antimicrobial iodine. After installation, check water flow and verify that the unit is operational. Create a check-off sheet to monitor volume of water used through unit. Shut off water to the dental unit when not in use.

The lifetime of this product is 240 liters or 365 calendar days (if water usage records are not monitored) or if iodine levels fall to less than 0.5ppm. After this period, or level is reached, immediately replace the spent Hu-Friedy Waterline Filter with a new or replacement unit. To install the new unit, follow the instructions above for installing a new product.

STORAGE & DISPOSAL

Storage: Store at room temperature in original packaging, away from children, pets and direct sunlight. Avoid excessive heat. Do not allow cartridge to freeze. Do not store in open or unlabeled containers.

Dispose: Securely wrap cartridge in newspaper and discard in trash.

NOTE: Adherence to CDC/ADA guidelines for aseptic procedures including a 2-minute morning flush, a 20-30 second flush between patients, and hand piece sterilization must be continued after the installation of this product.

Troubleshooting: If the water flow after product installation is slow or inhibited, check to ensure that the product was installed properly or if the product has been plugged by debris. For bottle installation ensure that the air supply is turned on at the maximum level (45 psi). In addition, ensure that the product is at the bottom of the water bottle. If there is an air leak from fittings, contact your dental supplier for replacement parts or service. Discoloration of water is caused by release. Initial flush of excess iodine. Operate water until flow appears clear.

Once the product is in place, the supply water will pass through the iodinated resin bed, which will continuously release low levels of antimicrobial iodine. After installation, check water flow and verify that the unit is operational. Create a check-off sheet to monitor volume of water used through unit. Shut off water to the dental unit when not in use.

Dispose: Securely wrap cartridge in newspaper and discard in trash.

NZ736/0618